



Since 1947

SEAL



ELECTROSURGICAL SOLUTIONS

Powered by Lamidey Noury Medical

Made in France

SEAL

The most versatile, efficient and affordable unit



Full electrocautery units 370W

Complete electrocautery unit with a wide range of monopolar and bipolar effects.



Vessel sealing system

Safe and efficient up to 7mm vessel with validation signal after sealing process.



PLASMA EDGE

The gold standard of bipolar plasma resection power dedicated for bipolar TUR and hystero-resection procedure.



Thermocision

A unique system which can Seal, Dissect and Cut in one instrument dedicated for Thyroid and delicate procedure.



Endosection

Specific waveform dedicated to endoscopic polypectomy.




Argon system

Argon Plasma Coagulation with the optional Argon module system.


Monopolar and Bipolar Accessories

Application:
↳ All Specialities

Reusable



Reusable




Bipolar/Monopolar



Single use



Single use



Cables



Patient Plate

Monopolar Pencil

Forceps/Cables

Endosection and Argon Capability

Application:
ENDOSCOPY
↳ Gastroenterology
OPEN SURGERY
↳ Liver coagulation



Argon capability with external argon module



Endosection power dedicated for polypectomy



Vessel sealing Forceps

Application:

- ↳ Visceral Surgery
- ↳ Gynecology
- ↳ Urology



Thermoclamp
16cm

Thermocut
5mm

Thermocut
10mm

Thermoclamp
22cm



Reusable
and autoclavable



Single use blade



Insulated jaws



Dissect, Seal and Cut




Application:

ENT

- ↳ Total Thyroidectomy
- ↳ Parotidectomy
- ↳ Carotid Access

BREAST

- ↳ Axillary Node Dissection



Reduced Thermal Damage

The thermal damage of a surgical instrument is caused by 2 critical elements:

- time applied on tissues
- energy delivered

The Thermocision technology is designed to reduce significantly both elements

Bipolar Plasma Resection System

Complete single use electrode range



UROLOGY LINE



Thin Loop



Thick Loop



Super Thick



Vaporization



Ball



Needle



Enucleation

GYNECOLOGY LINE



Slim loop



Cutting loop



Straight loop



Vaporization



Ball



Needle




Twizzle 5Fr

WORKING ELEMENT SOLUTION



Stand alone working elements compatible to usual resectoscopes

COMPLETE SOLUTION



Double flow sheaths quick lock design with dedicated working element



Since 1947

Designed and manufactured by Lamidey Noury



French manufacturer of electrosurgical solutions since 1947

Lamidey Noury Medical

Z.A. Les Godets 3 rue des Petits Ruisseaux - F 91370 Verrières le Buisson FRANCE

Tel. : +33 (0)1 69 20 69 69 - Fax : +33 (0)1 60 13 97 47

info@lamidey-noury.fr - export@lamidey-noury.fr

www.lamidey-noury.com



Made in France