



ARTHROSCOPY PUNCHES

POWERFUL PINLESS HINGE DESIGN

► Top Sellers



STELLA **LIGHT**
Coated light weighed aluminium/steel handle



STELLA **STANDARD**
Stainless steel handle



STELLA **LOCK**
Stainless steel handle with slide lock

ARTHROSCOPY PUNCHES



	Description	Tip profile in mm (A)	Bite width in mm (B)	Tip width in mm (C)	STELLA STANDARD	STELLA LIGHT
WIDE PUNCH LOW PROFILE						
	Straight medium	1.90	3.20	5.00	578-020-00	577-020-00
	Curved left 12°; medium	1.90	3.20	5.00	578-035-03	577-035-03
	Curved right 12°; medium	1.90	3.20	5.00	578-040-04	577-040-04
	Curved up 15°; medium	1.90	3.20	5.00	578-030-02	577-030-02
WIDE PUNCH LOW PROFILE; UPBITER						
	Upbiter straight; medium	1.90	3.20	5.00	578-025-10	577-025-10
WIDE PUNCH STANDARD PROFILE						
	Straight medium	2.30	3.20	5.00	578-120-00	577-120-00
	Curved up 15°; medium	2.30	3.20	5.00	578-130-02	577-130-02
WIDE PUNCH STANDARD PROFILE; UPBITER						
	Upbiter straight; medium	2.30	3.20	5.00	578-125-10	577-125-10
	Upbiter curved left 12°; medium	2.30	3.20	5.00	578-135-13	577-135-13
	Upbiter curved right 12°; medium	2.30	3.20	5.00	578-140-14	577-140-14
BASKET PUNCH SLIM						
	Straight	1.90	1.70	3.10	578-210-00	577-210-00
	Curved up 15°	1.90	1.70	3.10	578-215-10	577-215-10
BASKET PUNCH SLIM; UPBITER						
	Upbiter straight	1.90	1.70	3.10	578-210-10	577-210-10
	Upbiter curved left 12°	1.90	1.70	3.10	578-220-13	577-220-13
	Upbiter curved right 12°	1.90	1.70	3.10	578-225-14	577-225-14
BASKET PUNCH SCOOP STYLE						
	Straight; small	2.10	1.60	3.20	578-310-00	577-310-00
	Straight; small; wide basket	2.10	1.60	4.00	578-335-00	577-335-00
	Straight; medium	2.30	1.80	3.60	578-350-00	577-350-00
	Straight; large; narrow basket	2.60	2.20	3.60	578-355-00	577-355-00
	Curved left 12°; small	2.10	1.60	3.20	578-365-03	577-365-03
	Curved right 12°; small	2.10	1.60	3.20	578-370-04	577-370-04
	Curved up 15°; small	2.10	1.60	3.20	578-320-02	577-320-02
	Curved down 15°; small	2.10	1.60	3.20	578-360-05	577-360-05
BASKET PUNCH SCOOP STYLE; UPBITER						
	Upbiter straight; small	2.10	1.70	3.20	578-315-10	577-315-10
	Upbiter straight; medium	2.30	2.20	4.10	578-340-10	577-340-10
	Upbiter curved left 12°; small	2.10	1.70	3.20	578-325-13	577-325-13
	Upbiter curved right 12°; small	2.10	1.70	3.20	578-330-14	577-330-14
	Upbiter curved up 15°; medium	2.30	2.20	4.10	578-345-02	577-345-02

ARTHROSCOPY PUNCHES & GRASPERS



Description	Tip profile in mm (A)	Bite width in mm (B)	Tip width in mm (C)	STELLA STANDARD	STELLA LIGHT
BASKET PUNCH POSTERIOR					
Straight	2.30	2.20	4.00	578-385-00	577-385-00
Curved left 12°	2.30	2.20	4.00	578-385-03	577-385-03
Curved right 12°	2.30	2.20	4.00	578-385-04	577-385-04
Curved up 15°	2.30	2.20	4.00	578-385-02	577-385-02
BASKET PUNCH POSTERIOR; UPBITER					
Upbiter straight	2.30	2.20	4.00	578-385-10	577-385-10
OVAL PUNCH					
Straight	1.90	4.40	6.70	578-410-00	577-410-00
Curved left 12°	1.90	4.40	6.70	578-425-03	577-425-03
Curved right 12°	1.90	4.40	6.70	578-430-04	577-430-04
Curved up 15°	1.90	4.40	6.70	578-420-02	577-420-02
OVAL PUNCH; UPBITER					
Upbiter straight	1.90	4.40	6.70	578-415-10	577-415-10
SCISSORS PUNCH					
Straight	2.20	0.80	2.30	578-510-00	577-510-00
SCISSORS PUNCH; UPBITER					
Upbiter straight	2.20	0.80	2.30	578-515-10	577-515-10
Upbiter curved left 12°	2.20	0.80	2.30	578-520-13	577-520-13
Upbiter curved right 12°	2.20	0.80	2.30	578-525-14	577-525-14
BACKBITER PUNCH					
Straight; cutting left	3.80	2.40	5.60	578-555-30	577-555-30
Straight; cutting right	3.80	2.40	5.60	578-560-40	577-560-40

GRASPERS WITH SLIDE LOCK HANDLE

Description	Tip profile in mm	Tip width in mm	STELLA LOCK
GRASPER ALLIGATOR; MAXIMUM			
Straight	3.60	3.90	578-720-00
GRASPER ALLIGATOR; FINE			
Straight; fine	3.40	3.10	578-715-00
GRASPER; 1 x 2 TEETH			
Straight; fine	4.00	3.90	578-710-00
GRASPER; CUPPED			
Straight	3.30	2.80	578-730-00

For additional models refer to our Arthroscopy catalogue.



· WORLD WIDE DISTRIBUTION NETWORK ·



TONTARRA Medizintechnik GmbH

Daimlerstr. 15 · 78573 Wurmlingen/Germany

Tel: +49(0) 74 61/9 65 76 -0

Fax: +49(0) 74 61/9 65 76 -26

info@tontarra.de

www.tontarra.de

All rights are reserved. Technical alterations are possible. This folder may be used for no other purpose than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse, we retain the rights to recall our catalogs and pricelists and to take legal actions.
999-311-01E · 04.2019